


Welcome! Follow these steps to set up and use your Tyto device.


SET UP

You'll need a mobile phone or tablet that has internet access through Wi-Fi or a personal hotspot.

- 1 Download the TytoCare app from the App Store (iOS) or Google Play (Android) onto your smartphone or tablet.
- 2 Open the TytoCare app. Sign into your account with your email address and the temporary password you received from TytoCare. Be sure to change your password.
- 3 Turn on your Tyto device by pressing the button on the top of the device. Follow the instructions on your device to pair it with the TytoCare app.




LEARN

- 1 **Complete Tyto Academy**
Follow the Tyto tutorials to learn how to use your Tyto device. To start, tap the  icon in the upper right corner of the TytoCare app.
- 2 **Try It Out with Our Help**
Our support team will contact you about a week before your upcoming virtual visit. You'll practice using the Tyto device by taking and submitting a temperature exam.

PERFORM EXAM

TytoCare exams are used to share measurements with your provider in two different situations.

- **Before a Scheduled Virtual Visit**
Perform the exam one day before a scheduled visit.
- **For Non-Emergency Health Concerns**
The exams you submit are not monitored continuously. For follow-up, you **MUST** contact your clinic.



JOIN VISIT

To join a TytoCare virtual visit:

- 1 Log in to the TytoCare app and select **See a Clinician Online**. Then select your clinic.
- 2 When you are asked to join, select **Join** and enter the meeting.